

For Research Use Only

GSDMB Monoclonal antibody, PBS Only



Catalog Number: 67527-1-PBS

Basic Information

Catalog Number:

67527-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28332

GenBank Accession Number:

BC025682

GeneID (NCBI):

55876

UNIPROT ID:

Q8TAX9

Full Name:

gasdermin B

Calculated MW:

411 aa, 47 kDa

Observed MW:

45-50 kDa

Purification Method:

Protein G purification

CloneNo.:

2F12C4

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



MCF-7 cells were subjected to SDS PAGE followed by western blot with 67527-1-Ig (GSDMB antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67527-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed HepG2 cells using GSDMB antibody (67527-1-Ig, Clone: 2F12C4) at dilution of 1:800 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 67527-1-PBS in a different storage buffer formulation.