

For Research Use Only

STK32A Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67528-3-PBS

Basic Information

| | | |
|---|---|--|
| Catalog Number: 67528-3-PBS | GenBank Accession Number: BC021666 | Purification Method: Protein A Magarose purification |
| Size: 1 mg/ml | GeneID (NCBI): 202374 | CloneNo.: 2B6E5 |
| Source: Mouse | UNIPROT ID: Q8WU08 | |
| Isotype: IgG2a | Full Name: serine/threonine kinase 32A | |
| Immunogen Catalog Number: AG10942 | Calculated MW: 396aa,46 kDa; 166aa,20 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50267-2, STK32A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67528-1-PBS. Detection antibody: 67528-3-PBS. Standard: Ag10942. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP50267-3, STK32A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67528-2-PBS. Detection antibody: 67528-3-PBS. Standard: Ag10942. Range: 1.563-200 ng/mL.