

For Research Use Only

AMT Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67532-3-PBS



Basic Information

Catalog Number: 67532-3-PBS	GenBank Accession Number: BC007546	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 275	CloneNo.: 1F11E3
Source: Mouse	UNIPROT ID: P48728	
Isotype: IgG1	Full Name: aminomethyltransferase	
Immunogen Catalog Number: AG26644	Calculated MW: 44 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50403-1, AMT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67532-2-PBS. Detection antibody: 67532-3-PBS. Standard:Ag26644. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50403-2, AMT Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67532-3-PBS. Detection antibody: 67532-1-PBS. Standard:Ag26644. Range: 3.125-100 ng/mL.