

For Research Use Only

RALYL Monoclonal antibody, PBS Only



Catalog Number: 67534-1-PBS

Basic Information

Catalog Number: 67534-1-PBS	GenBank Accession Number: BC031090	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 138046	CloneNo.: 1F3F4
Source: Mouse	UNIPROT ID: Q86SE5	
Isotype: IgG1	Full Name: RALY RNA binding protein-like	
Immunogen Catalog Number: AG11676	Calculated MW: 291 aa, 32 kDa	
	Observed MW: 35-40 kDa	

Applications

Tested Applications:
WB, IF, ELISA

Species Specificity:
Human, mouse, rat, pig

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67534-1-Ig (RALYL antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67534-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using 67534-1-Ig (RALYL antibody), at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67534-1-PBS in a different storage buffer formulation.