

For Research Use Only

# TNFSF8 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67543-1-PBS



## Basic Information

<b>Catalog Number:</b> 67543-1-PBS	<b>GenBank Accession Number:</b> BC093630	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 944	<b>CloneNo.:</b> 2E10H5
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> P32971	
<b>Isotype:</b> IgG2a	<b>Full Name:</b> tumor necrosis factor (ligand) superfamily, member 8	
<b>Immunogen Catalog Number:</b> AG12726	<b>Calculated MW:</b> 234 aa, 26 kDa	

## Applications

**Tested Applications:**  
Cytometric bead array, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP50555-3, TNFSF8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67543-1-PBS. Detection antibody: 67543-4-PBS. Standard: Ag12726. Range: 6.25-100 ng/mL