

For Research Use Only

GRAMD1B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67551-4-PBS



Basic Information

Catalog Number:

67551-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG20670

GenBank Accession Number:

BC107480

GeneID (NCBI):

57476

UNIPROT ID:

Q3KR37

Full Name:

GRAM domain containing 1B

Calculated MW:

738 aa, 85 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D1A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50885-2, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-4-PBS. Detection antibody: 67551-3-PBS. Standard: Ag20670. Range: 0.049-3.125 ng/mL