

For Research Use Only

GRAMD1B Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67551-7-PBS



Basic Information

Catalog Number: 67551-7-PBS	GenBank Accession Number: BC107480	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 57476	CloneNo.: 1H9H2
Source: Mouse	UNIPROT ID: Q3KR37	
Isotype: IgG2a	Full Name: GRAM domain containing 1B	
Immunogen Catalog Number: AG20670	Calculated MW: 738 aa, 85 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50885-4, GRAMD1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67551-7-PBS. Detection antibody: 67551-3-PBS. Standard: Ag20670. Range: 1.563-50 ng/mL