

For Research Use Only

DDX58 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67556-3-PBS



Basic Information

Catalog Number: 67556-3-PBS	GenBank Accession Number: BC132786	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 23586	CloneNo.: 2F10G6
Source: Mouse	UNIPROT ID: O95786	
Isotype: IgG1	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 58	
Immunogen Catalog Number: AG18591	Calculated MW: 925 aa, 106 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50370-1, RIG-1/DDX58 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67556-2-PBS. Detection antibody: 67556-3-PBS. Standard: Ag18591. Range: 0.391-50 ng/mL