

For Research Use Only

GRP75 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67563-2-PBS



Basic Information

Catalog Number: 67563-2-PBS	GenBank Accession Number: BC000478	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3313	CloneNo.: 1D5G2
Source: Mouse	UNIPROT ID: P38646	
Isotype: IgG1	Full Name: heat shock 70kDa protein 9 (mortalin)	
Immunogen Catalog Number: AG7125	Calculated MW: 74 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50590-1, GRP75 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67563-2-PBS. Detection antibody: 67563-3-PBS. Standard: Ag7125. Range: 0.391-100 ng/mL.