

For Research Use Only

# RPS14 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67566-4-PBS



## Basic Information

Catalog Number:

67566-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30032

GenBank Accession Number:

BC003401

GeneID (NCBI):

6208

UNIPROT ID:

P62263

Full Name:

ribosomal protein S14

Calculated MW:

16 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2H7D3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50978-2, RPS14 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67566-4-PBS. Detection antibody: 67566-5-PBS. Standard:Ag30032. Range: 1.563-100 ng/mL.