

For Research Use Only

HSD3B2 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67572-2-PBS



Basic Information

Catalog Number:

67572-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG4291

GenBank Accession Number:

BC038419

GeneID (NCBI):

3284

UNIPROT ID:

P26439

Full Name:

hydroxy-delta-5-steroid
dehydrogenase, 3 beta- and steroid
delta-isomerase 2

Calculated MW:

42 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1F8B2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51043-1, HSD3B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67572-2-PBS. Detection antibody: 67572-3-PBS. Standard:Ag4291. Range: 3.125-200 ng/mL.