

For Research Use Only

HSD3B2 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67572-3-PBS



Basic Information

Catalog Number: 67572-3-PBS	GenBank Accession Number: BC038419	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 3284	CloneNo.: 2D6E2
Source: Mouse	UNIPROT ID: P26439	
Isotype: IgG1	Full Name: hydroxy-delta-5-steroid dehydrogenase, 3 beta- and steroid delta-isomerase 2	
Immunogen Catalog Number: AG4291	Calculated MW: 42 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51043-1, HSD3B2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67572-2-PBS. Detection antibody: 67572-3-PBS. Standard:Ag4291. Range: 3.125-200 ng/mL.