

For Research Use Only

# Vaspin Monoclonal antibody

Catalog Number: 67573-1-Ig



## Basic Information

Catalog Number:

67573-1-Ig

Size:

2200 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30136

GenBank Accession Number:

BC040857

GeneID (NCBI):

145264

UNIPROT ID:

Q8IW75

Full Name:

serpin peptidase inhibitor, clade A  
(alpha-1 antiproteinase, antitrypsin),  
member 12

Calculated MW:

414 aa, 47 kDa

Observed MW:

47 kDa

Purification Method:

Protein G purification

CloneNo.:

2A1B6

Recommended Dilutions:

WB 1:2000-1:10000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

Human

Positive Controls:

WB: human adipose tissue, human placenta tissue

## Background Information

### Storage

Storage:

Store at -20°C.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Human adipose tissue were subjected to SDS PAGE followed by western blot with 67573-1-Ig (Vaspin antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.