

For Research Use Only

ANGPTL4 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67577-3-PBS

Basic Information

Catalog Number:

67577-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30070

GenBank Accession Number:

BC023647

GeneID (NCBI):

51129

UNIPROT ID:

Q9BY76

Full Name:

angiopoietin-like 4

Calculated MW:

45 kDa

Purification Method:

Protein A purification

CloneNo.:

2D9B9

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP50037-1, ANGPTL4 Monoclonal Matched Antibody Pair - BSA and Azide free, measured by sandwich ELISA. Capture antibody: 67577-2-PBS. Detection antibody: 67577-3-PBS (HRP conjugated). Standard: Ag30070. Range: 0.313-20 ng/mL