

For Research Use Only

SSBP1 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67582-5-PBS



Basic Information

Catalog Number:

67582-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29976

GenBank Accession Number:

BC000895

GeneID (NCBI):

6742

UNIPROT ID:

Q04837

Full Name:

single-stranded DNA binding protein

1

Calculated MW:

148 aa, 17 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1D3E11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50750-2, SSBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67582-4-PBS. Detection antibody: 67582-5-PBS. Standard: Ag29976. Range: 0.391-100 ng/mL.



Cytometric bead array standard curve of MP50750-4, SSBP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67582-5-PBS. Detection antibody: 67582-7-PBS. Standard: Ag29976. Range: 0.098-12.5 ng/mL.