

For Research Use Only

CALU Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67585-3-PBS



Basic Information

Catalog Number: 67585-3-PBS	GenBank Accession Number: BC013383	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 813	CloneNo.: 2C1D2
Source: Mouse	UNIPROT ID: O43852	
Isotype: IgG2a	Full Name: calumenin	
Immunogen Catalog Number: AG12073	Calculated MW: 315 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50698-1, CALU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67585-2-PBS. Detection antibody: 67585-3-PBS. Standard: Ag12073. Range: 0.098-25 ng/mL