

For Research Use Only

# CALU Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67585-4-PBS



## Basic Information

Catalog Number:

67585-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG12073

GenBank Accession Number:

BC013383

GeneID (NCBI):

813

UNIPROT ID:

O43852

Full Name:

calumenin

Calculated MW:

315 aa, 37 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

2B2B8

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50698-3, CALU Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67585-2-PBS. Detection antibody: 67585-4-PBS. Standard: Ag12073. Range: 0.195-100 ng/mL.