

For Research Use Only

# Prx5 Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67599-4-PBS



## Basic Information

Catalog Number:

67599-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG12262

GenBank Accession Number:

BC110983

GeneID (NCBI):

25824

UNIPROT ID:

P30044

Full Name:

peroxiredoxin 5

Calculated MW:

214 aa, 22 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1C1G10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



Cytometric bead array standard curve of MP50572-2, Prx5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67599-4-PBS. Detection antibody: 67599-3-PBS. Standard:Ag12262. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50572-3, Prx5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67599-4-PBS. Detection antibody: 67599-5-PBS. Standard:Ag12262. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50572-4, Prx5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67599-2-PBS. Detection antibody: 67599-4-PBS. Standard:Ag12262. Range: 0.098-100 ng/mL