

For Research Use Only

SDHB Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67600-2-PBS



Basic Information

Catalog Number:

67600-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29868

GenBank Accession Number:

BC007840

GeneID (NCBI):

6390

UNIPROT ID:

P21912

Full Name:

succinate dehydrogenase complex,
subunit B, iron sulfur (lp)

Calculated MW:

32 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2C5A5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50553-1, SDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67600-2-PBS. Detection antibody: 67600-3-PBS. Standard:Ag29868. Range: 1.563-100 ng/mL.



Cytometric bead array standard curve of MP50553-3, SDHB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67600-4-PBS. Detection antibody: 67600-2-PBS. Standard:Ag29868. Range: 6.25-100 ng/mL.