

For Research Use Only

CD63 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67605-3-PBS



Basic Information

Catalog Number: 67605-3-PBS	GenBank Accession Number: BC002349	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 967	CloneNo.: 1B2B2
Source: Mouse	UNIPROT ID: P08962	
Isotype: IgG1	Full Name: CD63 molecule	
Immunogen Catalog Number: AG19690	Calculated MW: 26 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50558-1, CD63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67605-2-PBS. Detection antibody: 67605-3-PBS. Standard:Ag19690. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP50558-2, CD63 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67605-4-PBS. Detection antibody: 67605-3-PBS. Standard:Ag19690. Range: 1.563-100 ng/mL