

For Research Use Only

# CREB3L1,OASIS Monoclonal antibody, PBS Only (Capture)

Catalog Number:67617-2-PBS



## Basic Information

Catalog Number:

67617-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG28042

GenBank Accession Number:

BC015781

GeneID (NCBI):

90993

UNIPROT ID:

Q96BA8

Full Name:

cAMP responsive element binding  
protein 3-like 1

Calculated MW:

57 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

1G6A2

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP50829-1, CREB3L1, OASIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67617-2-PBS. Detection antibody: 67617-3-PBS. Standard: Ag28042. Range: 0.391-100 ng/mL



Cytometric bead array standard curve of MP50829-3, CREB3L1, OASIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67617-2-PBS. Detection antibody: 67617-5-PBS. Standard: Ag28042. Range: 1.563-100 ng/mL