

For Research Use Only

CREB3L1,OASIS Monoclonal antibody, PBS Only (Detector)

Catalog Number:67617-3-PBS



Basic Information

Catalog Number: 67617-3-PBS	GenBank Accession Number: BC015781	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 90993	CloneNo.: 2A12D1
Source: Mouse	UNIPROT ID: Q96BA8	
Isotype: IgG1	Full Name: cAMP responsive element binding protein 3-like 1	
Immunogen Catalog Number: AG28042	Calculated MW: 57 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50829-1, CREB3L1, OASIS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67617-2-PBS. Detection antibody: 67617-3-PBS. Standard: Ag28042. Range: 0.391-100 ng/mL.