

For Research Use Only

CKAP5 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67631-2-PBS

Basic Information

Catalog Number: 67631-2-PBS	GenBank Accession Number: BC120869	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9793	CloneNo.: 1E1D5
Source: Mouse	UNIPROT ID: Q14008	
Isotype: IgG2b	Full Name: cytoskeleton associated protein 5	
Immunogen Catalog Number: AG28770		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50224-1, CKAP5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67631-2-PBS. Detection antibody: 67631-1-PBS. Standard:Ag28770. Range: 0.391-100 ng/mL.