

For Research Use Only

MAT2B Monoclonal antibody, PBS Only



Catalog Number: 67643-1-PBS

Basic Information

Catalog Number:

67643-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG8649

GenBank Accession Number:

BC005218

GeneID (NCBI):

27430

UNIPROT ID:

Q9NZL9

Full Name:

methionine adenosyltransferase II,
beta

Calculated MW:

36 kDa, 38 kDa

Observed MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

2B7G8

Applications

Tested Applications:

WB, Indirect ELISA, IHC

Species Specificity:

Human, Mouse, Rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67643-1-ig (MAT2B antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated Alpha Tubulin Monoclonal antibody (HRP-66031) as loading control. This data was developed using the same antibody clone with 67643-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 67643-1-ig (MAT2B antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0). This data was developed using the same antibody clone with 67643-1-PBS in a different storage buffer formulation.