

For Research Use Only

CD226 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67653-3-PBS

Basic Information

Catalog Number: 67653-3-PBS	GenBank Accession Number: BC074787	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 10666	CloneNo.: 2A10H12
Source: Mouse	ENSEMBL Gene ID: ENSG00000150637	
Isotype: IgG1	UNIPROT ID: Q15762	
Immunogen Catalog Number: AG12403	Full Name: CD226 molecule	
	Calculated MW: 336 aa, 39 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50115-1, CD226 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67653-2-PBS. Detection antibody: 67653-3-PBS. Standard:Ag12403. Range: 0.391-100 ng/mL.