

For Research Use Only

PROCR Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67658-4-PBS



Basic Information

Catalog Number: 67658-4-PBS	GenBank Accession Number: BC014451	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10544	CloneNo.: 1C3A11
Source: Mouse	UNIPROT ID: Q9UNN8	
Isotype: IgG3	Full Name: protein C receptor, endothelial (EPCR)	
Immunogen Catalog Number: AG25628	Calculated MW: 238 aa, 27 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50576-2, PROCR Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67658-4-PBS. Detection antibody: 67658-5-PBS. Standard: Ag25628. Range: 1.563-200 ng/mL.