

For Research Use Only

ETFa Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67661-3-PBS



Basic Information

Catalog Number:

67661-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG30257

GenBank Accession Number:

BC015526

GeneID (NCBI):

2108

UNIPROT ID:

P13804

Full Name:

electron-transfer-flavoprotein, alpha polypeptide

Calculated MW:

30 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1F10C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50564-1, ETFA Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67661-2-PBS. Detection antibody: 67661-3-PBS. Standard: Ag30257. Range: 1.563-100 ng/mL.