

For Research Use Only

LRPPRC Monoclonal antibody, PBS Only (Capture/Detector)

Catalog Number: 67679-2-PBS



Basic Information

Catalog Number: 67679-2-PBS	GenBank Accession Number: BC050311	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 10128	CloneNo.: 3C8G4
Source: Mouse	UNIPROT ID: P42704	
Isotype: IgG1	Full Name: leucine-rich PPR-motif containing	
Immunogen Catalog Number: AG30261	Calculated MW: 1394 aa, 158 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50550-1, LRPPRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67679-2-PBS. Detection antibody: 67679-3-PBS. Standard: Ag30261. Range: 3.125-100 ng/mL.



Cytometric bead array standard curve of MP50550-2, LRPPRC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67679-4-PBS. Detection antibody: 67679-2-PBS. Standard: Ag30261. Range: 0.781-100 ng/mL.