

For Research Use Only

NOX4 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67681-4-PBS



Basic Information

Catalog Number: 67681-4-PBS	GenBank Accession Number: BC040105	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 50507	CloneNo.: 2C4E2
Source: Mouse	UNIPROT ID: Q9NPH5	
Isotype: IgG1	Full Name: NADPH oxidase 4	
Immunogen Catalog Number: AG6176	Calculated MW: 67 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51046-2, NOX4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67681-4-PBS. Detection antibody: 67681-5-PBS. Standard:Ag6176. Range: 0.391-25 ng/mL