

For Research Use Only

ALPK1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67689-3-PBS

Basic Information

Catalog Number: 67689-3-PBS	GenBank Accession Number: BC060780	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 80216	CloneNo.: 3A6H4
Source: Mouse	UNIPROT ID: Q96QP1	
Isotype: IgG1	Full Name: alpha-kinase 1	
Immunogen Catalog Number: AG5320	Calculated MW: 139 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50186-1, ALPK1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67689-3-PBS. Detection antibody: 67689-4-PBS. Standard: Ag5320. Range: 1.563-100 ng/mL.