

For Research Use Only

NDUFB8 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67690-2-PBS

Basic Information

Catalog Number: 67690-2-PBS	GenBank Accession Number: BC000466	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 4714	CloneNo.: 1H9G11
Source: Mouse	UNIPROT ID: O95169	
Isotype: IgG2a	Full Name: NADH dehydrogenase (ubiquinone) 1 beta subcomplex, 8, 19kDa	
Immunogen Catalog Number: AG6569	Calculated MW: 22 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50307-1, NDUFB8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67690-1-PBS. Detection antibody: 67690-2-PBS. Standard: Ag6569. Range: 0.098-6.25 ng/mL