

For Research Use Only

DLD Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67702-2-PBS



Basic Information

Catalog Number: 67702-2-PBS	GenBank Accession Number: BC018696	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1738	CloneNo.: 2A7A11
Source: Mouse	UNIPROT ID: P09622	
Isotype: IgG2b	Full Name: dihydrolipoamide dehydrogenase	
Immunogen Catalog Number: AG9680	Calculated MW: 509 aa, 54 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50544-1, DLD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67702-2-PBS. Detection antibody: 67702-3-PBS. Standard:Ag9680. Range: 3.125-100 ng/mL



Cytometric bead array standard curve of MP50544-2, DLD Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67702-2-PBS. Detection antibody: 67702-4-PBS. Standard:Ag9680. Range: 1.563-100 ng/mL