

For Research Use Only

RHAG Monoclonal antibody, PBS Only



Catalog Number: 67714-1-PBS

Basic Information

Catalog Number:

67714-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG30130

GenBank Accession Number:

BC012605

GeneID (NCBI):

6005

UNIPROT ID:

Q02094

Full Name:

Rh-associated glycoprotein

Calculated MW:

409 aa, 44 kDa

Observed MW:

50-60 kDa

Purification Method:

Protein G purification

CloneNo.:

1C1C4

Applications

Tested Applications:

WB, Indirect ELISA, IHC, IF

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




human placenta tissue were subjected to SDS PAGE followed by western blot with 67714-1-Ig (RHAG antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67714-1-PBS in a different storage buffer formulation.




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 67714-1-Ig (RHAG antibody) at dilution of 1:4000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9). This data was developed using the same antibody clone with 67714-1-PBS in a different storage buffer formulation.



Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using 67714-1-Ig (RHAG antibody) at dilution of 1:4000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer(pH9). This data was developed using the same antibody clone with 67714-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed human placenta tissue using RHAG antibody (67714-1-Ig, Clone: 1C1C4) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67714-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed human placenta tissue using RHAG antibody (67714-1-Ig, Clone: 1C1C4) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L). This data was developed using the same antibody clone with 67714-1-PBS in a different storage buffer formulation.