

For Research Use Only

CEP290 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67716-3-PBS

Basic Information

| | | |
|---|---|---|
| Catalog Number: 67716-3-PBS | GenBank Accession Number: BC008641 | Purification Method: Protein G purification |
| Size: 1 mg/ml | GeneID (NCBI): 80184 | CloneNo.: 2H4G5 |
| Source: Mouse | UNIPROT ID: O15078 | |
| Isotype: IgG1 | Full Name: centrosomal protein 290kDa | |
| Immunogen Catalog Number: AG17457 | Calculated MW: 2479 aa, 290 kDa | |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50156-1, CEP290 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67716-2-PBS. Detection antibody: 67716-3-PBS. Standard: Ag17457. Range: 1.563-100 ng/mL.