

For Research Use Only

# WWOX Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67725-2-PBS



## Basic Information

Catalog Number:

67725-2-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7969

GenBank Accession Number:

BC003184

GeneID (NCBI):

51741

UNIPROT ID:

Q9NZC7

Full Name:

WW domain containing  
oxidoreductase

Calculated MW:

47 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D11D7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51258-1, WWOX Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67725-2-PBS. Detection antibody: 67725-1-PBS. Standard:Ag7969. Range: 0.195-100 ng/mL.