

For Research Use Only

CTPS2 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67734-2-PBS

Basic Information

Catalog Number: 67734-2-PBS	GenBank Accession Number: BC034986	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 56474	CloneNo.: 1D10F9
Source: Mouse	UNIPROT ID: Q9NRF8	
Isotype: IgG1	Full Name: CTP synthase II	
Immunogen Catalog Number: AG30258	Calculated MW: 586 aa, 66 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP50097-1, CTPS2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67734-1-PBS. Detection antibody: 67734-2-PBS. Standard: Ag30258. Range: 0.195-50 ng/mL