

For Research Use Only

Mosquito AgBR1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67743-2-PBS

Basic Information

Catalog Number: 67743-2-PBS	GenBank Accession Number: XM_001844620	Purification Method: Protein G purification
Size: 1mg/ml	GeneID (NCBI): 6034353	CloneNo.: 2A1E3
Source: Mouse	UNIPROT ID: B0W7I8	
Isotype: IgG1	Full Name: AgBR1	
Immunogen Catalog Number: AG30168	Calculated MW: 49 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Mosquito

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50131-1, Mosquito AgBR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67743-2-PBS. Detection antibody: 67743-1-PBS. Standard: Ag30168. Range: 1.563-100 ng/mL