

For Research Use Only

KIF2C Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67752-3-PBS

Basic Information

Catalog Number: 67752-3-PBS	GenBank Accession Number: BC014924	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 11004	CloneNo.: 1C5B7
Source: Mouse	UNIPROT ID: Q99661	
Isotype: IgG1	Full Name: kinesin family member 2C	
Immunogen Catalog Number: AG29009	Calculated MW: 725 aa, 81 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50202-1, KIF2C Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67752-2-PBS. Detection antibody: 67752-3-PBS. Standard: Ag29009. Range: 0.781-100 ng/mL.