

For Research Use Only

ITCH Monoclonal antibody, PBS Only

Catalog Number: 67757-1-PBS



Basic Information

Catalog Number:

67757-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG14688

GenBank Accession Number:

BC011571

GeneID (NCBI):

83737

UNIPROT ID:

Q96J02

Full Name:

itchy E3 ubiquitin protein ligase homolog (mouse)

Calculated MW:

903 aa, 103 kDa

Observed MW:

103 kDa, 70 kDa

Purification Method:

Protein A purification

CloneNo.:

2E3G9

Applications

Tested Applications:

WB, Indirect ELISA, IF

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67757-1-Ig (ITCH antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 67757-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using ITCH antibody (67757-1-Ig, Clone: 2E3G9) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L), CL594-Phalloidin (red). This data was developed using the same antibody clone with 67757-1-PBS in a different storage buffer formulation.