

For Research Use Only

IGLC1 Monoclonal antibody, PBS Only (Capture)



Catalog Number: 67762-2-PBS

Basic Information

Catalog Number: 67762-2-PBS	GenBank Accession Number: BC007782	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 3537	CloneNo.: 3H2A7
Source: Mouse	Full Name: immunoglobulin lambda constant 1 (Mcg marker)	
Isotype: IgG2a	Calculated MW: 233 aa, 25 kDa	
Immunogen Catalog Number: AG14783		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50274-1, IGLC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67762-2-PBS. Detection antibody: 67762-3-PBS. Standard: Ag14783. Range: 3.125-50 ng/mL