

For Research Use Only

IGLC1 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67762-3-PBS

Basic Information

Catalog Number:

67762-3-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG14783

GenBank Accession Number:

BC007782

GeneID (NCBI):

3537

Full Name:

immunoglobulin lambda constant 1 (Mcg marker)

Calculated MW:

233 aa, 25 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

1A4A3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50274-1, IGLC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67762-2-PBS. Detection antibody: 67762-3-PBS. Standard: Ag14783. Range: 3.125-50 ng/mL



Cytometric bead array standard curve of MP50274-2, IGLC1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67762-4-PBS. Detection antibody: 67762-3-PBS. Standard: Ag14783. Range: 1.563-100 ng/mL