

For Research Use Only

SAR1B Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67765-5-PBS



Basic Information

Catalog Number:

67765-5-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG7718

GenBank Accession Number:

BC002847

GeneID (NCBI):

51128

UNIPROT ID:

Q9Y6B6

Full Name:

SAR1 homolog B (S. cerevisiae)

Calculated MW:

22 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

3B11C3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50803-2, SAR1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67765-4-PBS. Detection antibody: 67765-5-PBS. Standard:Ag7718. Range: 0.098-25 ng/mL



Cytometric bead array standard curve of MP50803-3, SAR1B Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67765-6-PBS. Detection antibody: 67765-5-PBS. Standard:Ag7718. Range: 0.781-100 ng/mL