

For Research Use Only

ALCAM Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67768-4-PBS



Basic Information

Catalog Number:

67768-4-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG16892

GenBank Accession Number:

BC137097

GeneID (NCBI):

214

UNIPROT ID:

Q13740

Full Name:

activated leukocyte cell adhesion molecule

Calculated MW:

583 aa, 65 kDa

Purification Method:

Protein A Magarose purification

CloneNo.:

3H2G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50574-2, ALCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67768-4-PBS. Detection antibody: 67768-5-PBS. Standard: Ag16892. Range: 6.25-200 ng/mL