

For Research Use Only

ALCAM Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67768-7-PBS



Basic Information

Catalog Number: 67768-7-PBS	GenBank Accession Number: BC137097	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 214	CloneNo.: 2F6D6
Source: Mouse	UNIPROT ID: Q13740	
Isotype: IgG1	Full Name: activated leukocyte cell adhesion molecule	
Immunogen Catalog Number: AG16892	Calculated MW: 583 aa, 65 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50574-4, ALCAM Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67768-5-PBS. Detection antibody: 67768-7-PBS. Standard: Ag16892. Range: 25-200 ng/mL