

For Research Use Only

GOLPH3 Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67777-2-PBS



Basic Information

Catalog Number:

67777-2-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG5443

GenBank Accession Number:

BC033725

GeneID (NCBI):

64083

UNIPROT ID:

Q9H4A6

Full Name:

golgi phosphoprotein 3 (coat-protein)

Calculated MW:

298 aa, 34 kDa

Purification Method:

Protein G Magarose purification

CloneNo.:

2D9E2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51157-1, GOLPH3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67777-2-PBS. Detection antibody: 67777-3-PBS. Standard: Ag5443. Range: 3.125-50 ng/mL