

For Research Use Only

GOLPH3 Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67777-3-PBS



Basic Information

Catalog Number: 67777-3-PBS	GenBank Accession Number: BC033725	Purification Method: Protein G Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 64083	CloneNo.: 3A10E7
Source: Mouse	UNIPROT ID: Q9H4A6	
Isotype: IgG1	Full Name: golgi phosphoprotein 3 (coat-protein)	
Immunogen Catalog Number: AG5443	Calculated MW: 298 aa, 34 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51157-1, GOLPH3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67777-2-PBS. Detection antibody: 67777-3-PBS. Standard: Ag5443. Range: 3.125-50 ng/mL