

For Research Use Only

RAB11B Monoclonal antibody, PBS Only



Catalog Number: 67780-1-PBS

Basic Information

Catalog Number:

67780-1-PBS

Size:

1 mg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG25429

GenBank Accession Number:

BC110081

GeneID (NCBI):

9230

UNIPROT ID:

Q15907

Full Name:

RAB11B, member RAS oncogene

family

Calculated MW:

24 kDa

Observed MW:

24 kDa

Purification Method:

Protein G purification

CloneNo.:

1G7D6

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

Human, mouse, rat

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 67780-1-Ig (RAB11B antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. The membrane was stripped and reblotted with HRP-conjugated GAPDH Monoclonal antibody (HRP-60004) as loading control. This data was developed using the same antibody clone with 67780-1-PBS in a different storage buffer formulation.