

For Research Use Only

CS Monoclonal antibody, PBS Only (Detector)

Catalog Number: 67784-5-PBS



Basic Information

Catalog Number: 67784-5-PBS	GenBank Accession Number: BC010106	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 1431	CloneNo.: 1C11B5
Source: Mouse	UNIPROT ID: O75390	
Isotype: IgG3	Full Name: citrate synthase	
Immunogen Catalog Number: AG9255	Calculated MW: 466 aa, 52 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50606-3, CS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67784-3-PBS. Detection antibody: 67784-5-PBS. Standard:Ag9255. Range: 0.098-100 ng/mL



Cytometric bead array standard curve of MP50606-4, CS Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67784-2-PBS. Detection antibody: 67784-5-PBS. Standard:Ag9255. Range: 0.098-100 ng/mL