

For Research Use Only

CD207 Monoclonal antibody, PBS Only (Detector)



Catalog Number: 67788-3-PBS

Basic Information

Catalog Number: 67788-3-PBS	GenBank Accession Number: BC022278	Purification Method: Protein A Magarose purification
Size: 1 mg/ml	GeneID (NCBI): 50489	CloneNo.: 2B5C5
Source: Mouse	UNIPROT ID: Q9UJ71	
Isotype: IgG2a	Full Name: CD207 molecule, langerin	
Immunogen Catalog Number: AG2479	Calculated MW: 328 aa, 37 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50310-1, CD207 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67788-2-PBS. Detection antibody: 67788-3-PBS. Standard:Ag2479. Range: 0.391-100 ng/mL.