

For Research Use Only

AIF Monoclonal antibody, PBS Only (Capture)

Catalog Number: 67791-2-PBS



Basic Information

Catalog Number: 67791-2-PBS	GenBank Accession Number: BC111065	Purification Method: Protein G purification
Size: 1 mg/ml	GeneID (NCBI): 9131	CloneNo.: 10F1F2
Source: Mouse	UNIPROT ID: O95831	
Isotype: IgG1	Full Name: apoptosis-inducing factor, mitochondrion-associated, 1	
Immunogen Catalog Number: AG12400	Calculated MW: 609 aa, 66 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP50508-1, AIF Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67791-2-PBS. Detection antibody: 67791-3-PBS. Standard: Ag12400. Range: 3.125-100 ng/mL.